



The Lindy and Norman Radow Sanctuary is equipped with a hearing loop system. If you wear a hearing aid or cochlear implants and it is equipped with a Telecoil (T-coil), you will be able to receive all audio which currently comes through the audio system, delivered directly through your hearing aid.

Often hearing aids have a switch that allows you to turn the T-coil on or some may automatically switch to T-coil when a signal is present. The switch position is usually indicated by a "T".

When you were fitted with your hearing aids, your audiologist should have told you if your hearing aid model has a T-coil. If you are not sure, it's best to call them and ask if they can confirm what hearing aid model you have. If it has a T-coil option, ask your audiologist how to switch to T-coil operation.



Congregation
ETZ CHAIM

creating meaningful Jewish experiences